

# UF-CPC 200



The Palas UF-CPC 200 is a butanol-based nanoparticle counter with high efficiency. It measures the number concentrations of ultrafine particles (UFP) in aerosols. Model 200 is designed for concentrations up to  $2 \cdot 10^6$  particles/cm<sup>3</sup>. This makes it especially suitable for low dilution measurement of high aerosol concentrations from particle generators, exhaust gas, or ambient air of exposed and heavily polluted areas.

In nephelometer mode, measurements up to  $2 \cdot 10^7$  particles/cm<sup>3</sup> are possible. The counter can be easily combined with the Palas size classifiers and thus directly allows the measurement of highly concentrated and size-selected aerosols, e.g., for efficiency measurement of HEPA filter media in test stands.

The patented evaporator and condensation module is maintenance-free. This allows continuous operating times of up to one year - unique.

## OPERATION PRINCIPLE

### CONDENSATION PARTICLE COUNTER FOR NUMBER CONCENTRATION UP TO $2 \cdot 10^6$ PARTICLES/CM<sup>3</sup>

The aerosol is fed directly to the evaporator unit by the internally adjustable and controlled diaphragm pump, saturated with the working liquid n-butanol. The butanol flows in a spiral trough inside the cylinder to the base of the evaporator. The unevaporated residue is pumped back into the reservoir by a second pump. This actively ensures permanent saturation of the evaporator while preventing deposits from forming on the inner walls.

In contrast to a control with a critical nozzle, contamination of the system cannot lead to a drop in volume flow. This is particularly important for long-term measurements. The user can also calibrate the volume flow.

The condensed particles are detected by an optical sensor, which determines the concentrations and size distribution of the condensed particles. This enables simple and efficient quality control.

The system is typically delivered with a cut-off D50 = 4.5 nm (measured with generated and selected NaCl particles). Others are adjustable and can be used optionally.

In the expert mode for research, the user can adjust various parameters to individual needs via the 7" touch screen.

The UF-CPC system supports a standardized interface with various protocol selection options, such as Modbus or ASCII protocol, for process monitoring applications. All measured data, including the corresponding settings, is stored directly on the device. Thus, the data can be accessed and visualized directly at any time.

## Extensions/Accessories

An air-conditioned weatherproof housing is available for the system.

## BENEFITS

- This novel and patented working fluid delivery system allows the user to quickly switch between butanol and isopropanol, for example, or even water.
- The UF-CPC 200 can count up to  $2 \cdot 10^6$  particles/cm<sup>3</sup> in single count mode, depending on the sensor in use (able to be switched out by the user).
- Integrated computer with 7" touch screen
- Intuitive user interface with sophisticated software for data analysis
- Integrated data logger
- Unlimited network compatibility that supports remote control and data storage on the Internet
- Integrated interface for process control applications

## DATASHEET

Measurement range (number $C_N$ )	$2 \cdot 10^6$ particles/cm <sup>3</sup> (single count mode)
Measurement range (size)	4 – 10,000 nm
Volume flow	0.9 l/min (butanol) 0.3 – 1 l/min (adjustable for research applications)
User interface	Touchscreen, 800 • 480 pixel, 7" (17.78 cm)
Data logger storage	4 GB
Software	PDAnalyze
Data acquisition	Digital, 20 MHz processor, 256 raw data channels
Light source	LED
Installation conditions	+10 – +30 °C (others on demand)
Accuracy	5 % (single count mode)
Response time	$t_{90} < 2.8$ s, $t_{90-10} < 2.0$ s
Operation liquid	1-butanol
Dimensions	290 • 240 • 350 mm (H • W • D)

## APPLICATIONS

- Aerosol research
- Testing of filters and air purifiers
- Environmental measurements
- Occupational exposure and workplace safety studies
- Inhalation and health effects studies
- Process monitoring
- Printer emission studies



Mehr Informationen:  
<https://www.palas.de/product/ufcpc200>