


# RESP-AER-METER SCIENTIFIC



Besides the detection of a potential superspreader, the Resp-Aer-Meter Scientific offers access to a wide range of other measurement data, such as number of particles, size and size distribution, time course, measurement times, and much more. In addition, many parameters, such as the measurement duration or limit values can be set individually.

## BENEFITS

- So-called "superemitters" can be identified in 30s thanks to a high number of particles
- Fast differentiation between infectious and less infectious Covid-19 carriers
- Measurement of the aerosol concentration and aerosol size in exhaled air
- Detection of particles from 145 nm to 10  $\mu\text{m}$
- Highest resolution, especially in the virus size range from approx. 145 nm to 1  $\mu\text{m}$
- Immediate evaluation and documentation of the measurement result

## APPLICATIONS

- Medical-scientific research, to answer questions such as
  - Do infectious respiratory diseases differ by exhaled particle size distribution?
  - Is an increased particle concentration in the breath due to a specific particle size range?
  - How do disease variations or personal characteristics (age, BMI, previous diseases...) influence the outcome?

## DATASHEET

Measuring principle	Optical light-scattering	Measurement (number $C_N$ )	range	0 – 20,000 particles/cm <sup>3</sup>
Measurement (size)	range	0.15 – 10 $\mu\text{m}$	Volume flow	9.5 l/min
User interface	Touchscreen, 800 • 480 pixel, 7" (17.78 cm)	Data acquisition		Digital, 20 MHz processor, 256 raw data channels
Power consumption	Approx. 200 W			